

A Type I lyophilised collagen absorbent primary dressing for Skin Donor sites

DonorDres® -An intact Type I collagen layer, backed by a hydrophilic sponge layer

- Intact Type I collagen layer facilitates healing^{1,2,3} and is hemostatic
- Provides immediate pain relief
- Collagen layer is biodegradable
- Hydrophilic sponge layer absorbs excess exudate

Indications:

- Skin donor sites, wide range of additional woundcare applications

Contraindications:

- Patients with known hypersensitivity to collagen
- The product is not sterile when package is damaged during transit
- DonorDres® is a sterile product which should not be resterilized

Ordering Information

CODE	DESCRIPTION	UNITS
DNDS1001-05	DonorDres® absorbent collagen dressing , 6cm x 7cm	Box of 5
DNDS1002-05	DonorDres® absorbent collagen dressing , 10cm x 10cm	Box of 5
DNDS1003-05	DonorDres® absorbent collagen dressing , 10cm x 15cm	Box of 5
DNDS1004-05	DonorDres® absorbent collagen dressing , 10cm x 20cm	Box of 5
DNDS1005-05	DonorDres® absorbent collagen dressing , 10cm x 30cm	Box of 5
DNDS1006-05	DonorDres® absorbent collagen dressing , 15cm x 30cm	Box of 5


Derived from Piscean sources

- Biocompatible as per EN ISO 10993 Standards
- Non-toxic, non-allergenic, non-immunogenic, non-pyrogenic

Storage:

Store in a dry place at temperatures not exceeding 25°C. Do not freeze.

Sterilised by irradiation with a 3 year shelf life.



References

¹ Hampton S, Understanding the pH balance in wound healing, JCN 2008;22:05

² Ponten B, Nordgaard JO. The Use of Collagen Film as a Dressing for Donor Areas in split skin grafting. Scand J Plast Reconstr Surg 1976; 10 : 237-40.

³ Krieg T. Collagen in the Healing Wound. WOUNDS. 1995;7(Suppl A):5A-12A



Manufactured by:

EUCARE PHARMACEUTICALS PRIVATE LIMITED

Plot No. AC 25-B, SIDCO Industrial Estate, Thirumudivakkam, Chennai - 600 044, India.
Tel: +91-44-4598 9955 / 9909 / 9904 / 9908 / 9903 / 9915, Fax: +91-44-2478 2516 / 4598 9914
E-mail: eucare@vsnl.com, Website : www.eucare.in / www.eucareindia.com

We create and surpass our own standards

DonorDres® - Registered Trade Mark in India.

EUCARE
the collagen people